

Certificate of Analysis




Website: <https://peptide.wumolab.com/>

Contact at peptide@wumolab.com

Product Name	BPC-157 5 mg/vial	Received Date	June 18, 2026
Batch No.	WUBC0616-001	Reported Date	June 25, 2026
Brand	WUMO Peptide	Manufacturer	SUZHOU XILIN TECHNOLOGY CO.,LTD
Sample description	White lyophilized powder in labeled vial	Tests requested	Identity and purity assessment of peptide vial sample


Specification	Result	Method	Status
Identity	BPC-157	ESI-MS	✓ Conforms
Fill Amount / Label Claim	5 mg/vial	Filling Record	✓ Conforms
Purity	99.31%	RP-HPLC (220 nm)	✓ Conforms

HPLC Chromatogram (Retention Time: 8.262 min)



No.	Peak Name	Retention Time min	Area mAU*min	Height mAU	Relative Area %	Relative Height %	Amount
n.a.	Component I	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
1		7.898	1.475	2.789	0.67	1.69	n.a.
2		8.262	157.577	412.549	99.31	96.97	n.a.
3		10.277	1.681	0.083	0.02	1.34	n.a.
Total			160.733	416.211	100	100.00	

Mass Spectrum (ESI-MS)



Remarks: For laboratory research use only. Not for human or veterinary use.

Checked By

Andy Gong

This report summarizes the test results for the submitted sample. Purity was analyzed by RP-HPLC, and identity was assessed by ESI-MS. The results apply only to the tested batch/sample and are provided for research and quality-control reference only.