

Certificate of Analysis




Website: <https://peptide.wumolab.com/>

Contact at peptide@wumolab.com


Product Name	GHK-Cu 50 mg/vial	Received Date	June 18, 2026
Batch No.	WUCU0616-001	Reported Date	June 25, 2026
Sample description	Blue lyophilized powder in labeled vial	Tests requested	Identity and purity assessment of peptide vial sample

Specification	Result	Method	Status
Identity	GHK-Cu	ESI-MS	✓ Conforms
Fill Amount / Label Claim	50 mg/vial	Filling Record	✓ Conforms
Purity	99.969%	RP-HPLC (220 nm)	✓ Conforms

HPLC Chromatogram (Retention Time: 4.501 min)



Mass Spectrum (ESI-MS)



Remarks: For laboratory research use only. Not for human or veterinary use.

Checked By

Andy Gong

This report summarizes the test results for the submitted sample. Purity was analyzed by RP-HPLC, and identity was assessed by ESI-MS. The results apply only to the tested batch/sample and are provided for research and quality-control reference only.